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**Instruction manual EN**

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# ***Sensorturn*** pro

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with the enclosed description fulfills the requirements of the Low Voltage Directive 2006/95/EC  
(valid to 19 April 2016) the Low Voltage Directive 2014/35/EU (valid from 20 April 2016),  
the regulations of the EMC Directive 2014/30/EC and the eco design Directive 2009/125/EC.

Hereby, we certify that our products, regardless of the production location, RoHS compliant and  
fulfill the directive 2011/65/EC (revised version: directive 2002/95/EC).

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

- 1 - Knob, Off - On and Rotation Speed
- 2 - Knob, Constant - Time
- 3 - IR-Sensor
- 4 - LED Constant
- 5 - LED Time
- 6 - LED Footpedal
- 7 - Turntable
- 7a - Silicone Covering
- 7b - Fixing screw

- 8 - Power connector
  - 9 - Connector for foot pedal
  - 10 - Setting screw for IR-Sensor detection range
  - 11 - Multiswitch
- Switch 1: Enable foot pedal  
 Switch 2: Timer selection  
 (1 - 25 sec. or 5 - 125 sec.)  
 Switch 3: Enable foot pedal operation "Standard"

Read these instructions carefully to familiarize yourself with the product. Please retain these operating instruction for future reference.

**Use:**

The **Sensorturn pro** turntable has been developed for rational and sterile plating / inoculation of petri dishes. For operation without fatigue, the **Sensorturn pro** is equipped with an IR-Sensor and allows switching-on and -off of the turntable with movements of the hand. The **Sensorturn pro** is highly efficient, UV-resistant and flame-sterilizable.

EN

**SAFETY PRECAUTIONS:**

- **On unpacking the unit, check for possible transportation damages. If there is any visible damage, do not try to operate the device. (Contact our service representative immediately.)**
- **For safe operation the unit should be on a level surface. The ambient temperature should be +5° C to +50° C.**
- **The default supply must be connected to a voltage source of 100 - 240 V / 50/60 Hz. The device may only be connected to the mains power supply if the values for the operating voltage given above are identical to the available mains voltage.**
- **Protect against moisture.**
- **On ending operation, switch off the device with knob (1), (Off).**

## The range:

### Sensorturn pro

Art-No.: 7.002.000

includes IR-Sensor and connector for foot pedal  
14 - 210 rpm

5 operating modes:

IR-Sensor: Time / Constant

Foot pedal: Standard / Start-Stop / Time


Turn-table with silicone covering

Silicone mat for centring ring (14)

Screwdriver for IR-Sensor detection range (15)

Power connection (16)

Instruction manual and 2-year warranty

## 1. Operating location:

For safe operation the **Sensorturn pro** should be on a level surface. The ambient temperature should be +5° C to + 50° C.

## 2. Electrical Connection:

Insert the power connection (15) into the socket (8) on the back panel of the unit. The default supply must be connected to a voltage source of 100 - 240 V / 50/60 Hz. The device may only be connected to the mains power supply if the values for the operating voltage given above are identical to the available mains voltage.

## 3. Switching on / off and rotation speed adjustment:

Switch the device on with the knob (1) and concurrently select the desired rotation speed (RPM). One or more LED's (4,5,6) light up depending on the selected operating mode. Switch the unit off by turning the knob (1) counter-clockwise to the Off-position.

### 3.1 LED's for operating modes:

LED Constant (4) lights up:

IR-Sensor operating mode "Constant"

LED Time (5) lights up:

IR-Sensor operating mode "Time"

LED Constant (4) and

LED Footpedal (6) light up:

Foot pedal operating mode "Start-Stop"

LED Time (4) and LED Footpedal (6) light up:

Foot pedal operating mode "Time"

LED Footpedal (6) lights:

Foot pedal operating mode "Standard"

Factory Setting: LED Constant (4) for IR-Sensor operating mode "Constant" lights up.

### 3.2 IR-Sensor operating mode "Constant"

Set the knob (2) to „Constant“ and ensure that the Multiswitch (11-1) "Enable foot pedal" at the rear side of the unit is in Off-position. To start the rotation, activate the IR-Sensor (3) by moving your hand past it at a distance of 5 to 50 mm.

The rotary movement of the turntable (7) is stopped by a second hand movement.

### 3.3 IR-Sensor operating mode "Time"

Set the knob (2) to „Time“ and ensure that the Multiswitch (11-1) "Enable foot pedal" at the rear side of the unit is in Off-position. To start the rotation, activate the IR-Sensor (3) by moving your hand past it at a distance of 5 to 50 mm. The turntable (7) stops automatically when the time set on the rotary switch (2) has run out.

### 3.4 Foot pedal operating mode “Start-Stop“:

Insert the cable of the foot pedal in the connector (9). Set the Multiswitch (11-1) “Enable foot pedal“ at the rear side of the unit to On-position and turn the knob (2) to “Constant“.

To start the rotation, press the foot pedal. The turntable (7) stops when the foot pedal is pressed for a second time.

Factory setting: Multiswitch 11-1 in Off-position

### 3.5 Foot pedal operating mode “Time“:

Insert the cable of the foot pedal in the connector (9). Set the Multiswitch (11-1) “Enable foot pedal“ at the rear side of the unit to On-position and turn the knob (2) to “Time“.

To start the rotation, press the foot pedal. The turntable (7) stops automatically when the time set on the rotary switch (2) has run out.

Factory setting: Multiswitch 11-1 in Off-position

### 3.6 Foot pedal operating mode “Standard“:

Insert the cable of the foot pedal in the connector (9). Set the Multiswitch (11-1) “Enable foot pedal“ and Multiswitch (11-3) “Enable foot pedal operation Standard“ at the rear side of the unit to On-position. The turntable (7) starts when the foot pedal is pressed. The foot pedal remains depressed for the duration of use. The turntable (7) stops upon release of the pedal.

Factory setting: Multiswitch 11-1 and 11-3 in Off-position

## 4. Detection range of the IR-Sensor:

By turning the setting screw (10) at the rear side of the unit (use the included screwdriver (15)), the detection range of the IR-Sensor can be adapted to special requirements. The detection range is shortened by turning it to the left; extended by turning it to the right.

Attention: It is possible to adjust the sensor range to 0 mm or more than 50 mm. Then the IR-Sensor is out of range and cannot be activated.

## 5. Choosing the time range for the operating mode “Time“:

When working in “Time“ operation the adjustable time range at knob (2) can be chosen between 1 - 25 sec. or 5 - 125 sec.

The time range can be chosen at Multiswitch (11-2) “Timer selection“ at the rear side of the unit.


Multiswitch 11-2 in Off-position: Time range 1 - 25 sec.

Multiswitch 11-2 in On-position: Time range 5 - 125 sec.

Factory setting: Multiswitch 11-2 in Off-position

## 6. Centering ring:

The centering ring can be used to center an individual type of petri dish on the turntable (7). To make the centering ring suitable for the needed petri dish size put the silicone mat for the centering ring (14) on the template printed in this instruction manual. Then put an empty petri dish upside down on the silicone mat and adjust it concentric with the help of the crosslines. Use a ball pen to mark the outside shape of the petri dish on the silicone

centering ring. Finally cut out the centering ring according to the marked circle and place the ring on the turntable (7).

### 7. Turntable cleaning and disassembly:


Check the unit is disconnected. The turntable (7) and silicone covering (7a) can be cleaned with customary commercial disinfectants. Additionally, it is possible to remove the turntable from the device to clean it separately. To remove the turntable proceed as follows: Remove the silicone covering (7a). Hold the turntable (7) with your hand and unscrew the fixing screw (7b) completely with a crosstip screwdriver. Turn approx. 4-5 revolutions in counter-clockwise direction. Now remove the turntable (7) from the device.

The disassembled turntable and the silicone covering can also be sterilized in an autoclave.

Note: If the turntable is sterilized in an autoclave the turntable may change it's colour slightly.

Reinstallation: Place the turntable (7) on the device, hold the turntable (7) with your hand and tight the fixing screw in clockwise direction.

Put back the silicone covering (7a) on the turntable (7).

### 8. Changing the turntable:

The included turntable for petri dishes up to a diameter of 100 mm can be interchanged with an larger turntable for petri dishes up to a diameter of 150 mm (optional, **Art.-No.: 7.000.150**).

To interchange the turntable proceed as follows: Remove the silicone covering (7a). Hold the turntable (7) with your hand and unscrew the fixing screw (7b) completely with a crosstip screwdriver. Turn approx. 4-5 revolutions in counter-clockwise direction. Now remove the turntable (7) from the device. Then, place the larger turntable on the device, hold the turntable with your hand and tight the fixing screw (7a) in clockwise direction. Put the larger silicone covering (included with the larger turntable **Art.-No.: 7.000.150**) on the turntable.

### 9. Warranty:

All **Sensorturn pro** are covered under our two-year manufacturer warranty against any manufacture defects in material and workmanship. The WLD-TEC warranty guarantees all **Sensorturn pro** under normal usage conditions and does not cover any damages as a direct result of user misuse or/and abuse. The warranty is void upon any unauthorized servicing, disassembly or modifications.

## Accessories:

### **Glass spatula**

in triangle form for inoculation

**Art.-No.: 6.000.370**

### **Stainless steel spatula**

in triangle form for inoculation

**Art.-No.: 6.000.371**

### **Tray**

mobile tray for instruments

**Art.-No.: 6.000.320**

### **Turntable Sensorturn pro, extra wide**

for petri dishes up to 150 mm diameter,  
with silicone covering

**Art.-No.: 7.000.150**

### **Tray, screw-on**

for right or left side, 3 places

**Art.-No.: 8.000.340**

### **Foot pedal made of stainless steel**

cable length 3 m

**Art.-No.: 6.000.402**

### **Foot pedal Mini, made of plastic**

cable length 3 m

**Art.-No.: 6.000.403**

### **Wireless radio frequency foot pedal**

(EU-countries only)

**Art.-No.: 7.000.404-RF**

## Troubleshooting guide

### **The green LED's (4,5,6) do not light up**

Check for correct connection and specification of the power adapter.  
Ensure that the original power adapter is used.

Specifications: 9V / DC, 1,3 A

Polarity: + —  -

### **One of the green LED's (4,5,6) is blinking fast**

This signal indicates a malfunction if the interior temperature has exceeded 70 °C.  
At a normal room temperature the unit is suited for continuous operation.  
In case of overtemperature. Allow sufficient time for the unit to cool down.

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
Firmenstempel / Company stamp

Armin Wegener  
Vice President Research & Development

## Technical data:

Technology Microprocessor

## Programs

IR-Sensor: Constant  
Time with timer (1 - 25 sec. or 5 - 125 sec., adjustable)

Foot pedal: Start-Stop  
Standard (rotating during pressed foot pedal)  
Time with timer (1 - 25 sec. or 5 - 125 sec., adjustable)

## Mechanical

Petri dish size: up to Ø 100 mm (Ø 150 mm optional)

Rotational direction: clockwise

Casing and operating controls: stainless steel / glass, UV and solvent resistant

Measurements (B x H x T): 109 x 59 x 130 mm

Weight: 860 g

Max. load: 400 g

## Electrical

Rotary speed range: 14 - 210 rpm

Detection range IR-Sensor: 5 - 50 mm, adjustable


Timer range "Timer" operation: 1 - 25 sec. or 5 - 125 sec., adjustable

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