**SteriMax** *basic*

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Instruction manual EN

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*Electrical equipment for measurement, control and laboratory use, EMC requirements*

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**Generic Emission Standard:**

*Electrical Equipment, class B, Group 1*

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*Safety requirements for electrical equipment for measurement, control, and laboratory use. Part 1: General requirements*

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*Safety requirements for electrical equipment for measurement, control, and laboratory use. Part 2-010: Particular requirements for laboratory equipment for the heating of materials*

### 3. Schutz vor optischer Strahlung (Risikofreie Gruppe) / *Security of optical radiation (Risk free group)*

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
*Photobiological safety of lamps and lamp systems*

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
**Notizen:**

Read these operating instructions carefully to familiarize yourself with the product before using. Keep the manual in a safe place for future reference.

Fig. A

- 1a - Start symbol - Touch panel
- 1b - Stop symbol - Touch panel
- 2 - Status display: sterilization progress
- 3 - Residual heat display
- 4 - Quartz glass annealing tube / annealing tube
- 4a - Sterilization orifice
- 5 - Warning symbol "Hot Surface"
- 6 - Magnets (inside) for optional tray

- 7 - Ventilation grilles
- 8 - Fan
- 9 - Drill hole for replacing annealing tube
- 10 - Connector for foot pedal
- 11 - Fuse holder
- 12 - Power socket
- 13 - Power switch
- 14 - Connector
- 15 - Power plug
- 16 - Power cable

Fig. B


Fig. C

## Table of contents

<b>The range</b> .....	18
<b>Use</b> .....	18
<b>Safety Precautions</b> .....	18
<b>1. Initial operation</b> .....	20
1.1 Quartz glass annealing tube .....	20
1.2 Power connection .....	20
<b>2. Operation</b> .....	20
2.1 Switching On / Off the device .....	20
2.2 Sterilization at start-up .....	20
<b>3. Sterilizing the instrument</b> .....	20
3.1 Starting the sterilization .....	21
3.2 Foot pedal connection and operation .....	21
3.3 Sterilization process .....	21
3.4 Cancellation of Sterilization .....	21
3.5 Thermocontrol .....	21
<b>4. Stand-by function</b> .....	22
<b>5. Residual heat display</b> .....	22
<b>6. Cleaning</b> .....	22
6.1 Cleaning and replacing the annealing tube .....	22
6.2 Cleaning the glass front .....	23
6.3 Cleaning the device .....	23
<b>7. Fan</b> .....	23
<b>8. Error display</b> .....	24
8.1 Over-temperature protection .....	24
8.2 IR-halogen lamp fault .....	24
8.3 Temperature sensor fault .....	24
<b>9. Warranty</b> .....	24
<b>Troubleshooting guide</b> .....	25
<b>Technical Data</b> .....	26
<b>Declaration of Conformity</b> .....	27
<b>Optional accessories</b> .....	28
<b>Spare parts</b> .....	28

## The range


**SteriMax basic - Inoculation Loop Sterilizer:**

**230V / AC Art. No.: 4.101.000**

**120V / AC Art. No.: 4.102.000**

**100V / AC Art. No.: 4.103.000**

- 1 quartz glass annealing tube (installed in the device)
- 1 inoculation loop holder (**R1**)
- 1 inoculation loop Ø 3 mm (**R1a**)
- Power cable (**R2**)
- 2 year warranty

## Use

The SteriMax basic infrared inoculating loop sterilizer is ideal for laboratory use for sterilizing inoculating loops, tweezers and micro-instruments made of heat resistant materials such as stainless steel or platinum-iridium alloys.

## Safety Precautions

- All users who have been assigned to use this device must read and understand this manual or have been instructed by a competent person in such a way that they can use this device without causing danger.
- When unpacking, check for possible shipping damage and do not operate if there is visible damage.
- Operate the unit only on a stable, non-slip and level surface.
- Do not operate the device near flammable liquids, materials or in potentially explosive areas.
- Never insert instruments or objects where flammable liquids or material adhere to these instruments or objects into the annealing tube (4, 4a).
- Operating conditions: Operating temperature 5°C to 40°C. Maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C.

- Make sure that the unit does not come in contact with water or other liquids as long as the plug is in the power socket, especially if you use it in the laboratory near water or liquids.
- Using a two-pin AC adapter / connecting cable is NOT recommended.

- Make sure that during operation the power cable (16) never gets wet or damp. Locate the cable so that it is not trapped, jammed or otherwise damaged. If the power cable (16) or the plug (15) is damaged, the cable must first be replaced before using this device.
- Do not cover the ventilation grilles (7) and cooling fan (8) on the device.
- Never insert plastic inoculating loops or objects with a temperature-resistance of less than 1200°C into the annealing tube (4, 4a).

- Never hold fingers or other parts of the body on or in the sterilization orifice (4a). Danger of burns!

- If the IR-halogen lamp is illuminated, do not look directly into the sterilization orifice (4a). Distance eye > 200 mm and distance skin > 30 mm. Never operate the device without the annealing tube in place (4)!

- Only use original quartz glass annealing tube.
- Only remove the annealing tube (4) after it has cooled down. Even some time after use, the sterilization orifice (4a) and the annealing tube (4) are still hot. Danger of burns! Allow the unit to cool (apt. 10 min.) and shut it down before cleaning, disinfection, maintenance or transport.
- After cleaning, let the annealing tube (4) dry first before reinstalling.
- Do not handle the annealing tube with bare hands. Grease or other dirt can be removed with a soft cloth with alcohol.
- After using the device or if it is unused for longer periods, turn off the power at the power switch (13).

- The upper part of the glass front can become hot if the device is in continuous use.

- In the event of faults or before opening the device, remove the mains plug (15) from the power socket.

- Keep all equipment and items that can be damaged by a magnetic field (such as credit & debit cards, data disks, mechanical watches, etc.) well away from the magnetic surfaces of the glass front (6).

For example the minimum distance to the glass front (6):  
credit & debit cards 32 mm, heart pacemaker 40 mm.

## 1. Initial operation

Please handle this sterilizer always with care as described in the following points of this manual.

### 1.1 Quartz glass annealing tube

Upon delivery, the annealing tube **(4)** is already installed in the sterilization orifice **(4a)**.

**Remove the protective foil and the adhesive warning label before first use. Never operate the device without the annealing tube (4) in place!**

### 1.2 Power connection

Plug the power cable **(16)** using the connector **(14)** into the power socket **(12)** on the back of the device. Insert the power plug **(15)** into a live power outlet suitable for your device (230V AC, 120V AC or 100V AC).

## 2. Operation

This section describes switching on and off the device and the sterilization at start-up.

### 2.1 Switching On / Off the device

1. Switch on the device using the power switch **(13)** on the rear panel **(switch position "I")**.
2. Switch off the device using the power switch **(13)** on the rear panel **(switch position "0")**.

Depending on the residual heat, after switching off, the fan **(8)** runs until the excess heat has been dissipated.

### 2.2 Sterilization at start-up

The sterilization at start-up starts automatically a sterilization after the device is switched on at the power switch **(13)**.

The blue start symbol **(1a)** lights up the first time, after completion of the sterilization at start-up.

## 3. Sterilizing the instrument

This section describes the touch operation of the glass front and the operation of the device for sterilizing inoculation loops or small instruments.

### 3.1 Starting the sterilization

Inserting an instrument (such as an inoculating loop with holder) into the orifice of the annealing tube **(4, 4a)**. By briefly touching the start symbol **(1a)** at the touch panel, the user can start the sterilization.

### 3.2 Foot pedal connection and operation

Alternatively, the sterilization process can also be started with a foot pedal, that is available optionally.

Insert the connection cable of the foot pedal into the socket **(10)** at the back of the unit. The sterilization progress can be started by operation of the foot pedal.

<b>Stainless steel foot pedal:</b>	<b>Art.-No.: 6.000.402</b>
<b>Foot pedal mini / plastic:</b>	<b>Art.-No.: 6.000.403</b>
<b>Wireless stainless steel foot pedal:</b>	<b>Art.-No.: 8.000.404-RF</b>

(Wireless for EU countries only)

### 3.3 Sterilization process

Hold the instrument in the annealing tube **(4)** until the sterilization process is finished (10 seconds). The status display **(2)** indicates the progress of sterilization due to the more increasing orange lighting of it. After 10 seconds, the sterilization will stop automatically and the status display fades out.

If the device is already heated by previous use, thermocontrol uses the residual heat and switching on the IR-halogen lamp will be delayed by up to 5 seconds.

Regardless of the unit temperature, the sterilizer is ready to use again after 1 - 3 seconds. The start symbol **(1a)** indicates this by lighting up blue.

Never hold fingers or other parts of the body on or in the sterilization orifice **(4a)**. Even some time after use, the sterilization orifice **(4a)** and the annealing tube **(4)** are still hot. Danger of burns!

If the IR-halogen lamp is illuminated, do not look directly into the sterilization orifice **(4a)**.

### 3.4 Cancellation of Sterilization

During the sterilization process, the stop symbol **(1b)** lights up. By lightly touching the symbol (analogous to starting the sterilization), the sterilization can be stopped. By operating the foot pedal (optional, **see section 3.2**) the sterilization can be stopped prematurely by pressing the foot pedal again.

### 3.5 Thermocontrol

Depending on the frequency of sterilization, the thermocontrol makes use of the residual heat, thus shortening the activation of the IR-halogen lamp automatically and dynamically by up to 5 seconds.

#### 4. Stand-by function

The stand-by time is 15 minutes.

While the stand-by time, the device switches off automatically after 15 minutes without operation to prevent the unwanted activation of a sterilization.

The fan **(8)** and the residual heat display **(3)** are active in standby mode until the excess heat has been dissipated.

#### 5. Residual heat display

The residual heat display **(3)** lights up red and warns that the sterilizing orifice is hot. The symbol is still displayed in the stand-by mode until the annealing tube has cooled down.

If the power supply of the unit is switched off at the power switch **(13)**, or the power cable is removed, the safety symbol cannot be displayed even if the annealing tube is still hot.

The fan continues to run for some time until the device has cooled down.

#### 6. Cleaning

**Before cleaning the quartz glass annealing tube, the device or the glass front:**

1. Let the appliance cool down until the residual heat display **(3)** does not more lights up.
2. Switch off the device using the power switch **(13)** on the rear panel **(switch position "0")**.

##### 6.1 Cleaning and replacing the annealing tube

The inside of the annealing tube should be cleaned regularly. Use the handle of the inoculating loop holder **(R1)** (included in the delivery) to push through the drill hole **(9)** on the rear of the device  and push the annealing tube forwards out through the sterilization orifice **(4a)**.

Use non-alkaline detergent for cleaning. Fill the annealing tube with appropriate cleaning solution and allow to stand for several hours to loosen burned-on residue.

(If necessary acidify in 5% hydrofluoric acid (HF) for 2-3 minutes). After cleaning the inside of the annealing tube, rinse well with distilled (deionized) water and allow to dry.

To lengthen the annealing tube's life, as a final step we recommend rubbing the outside with alcohol and then only touching it with gloves to re-insert it.

To install, slide the annealing tube with the closed end first into the sterilization orifice **(4a)** until it is fully inserted and is not protruding at the front. When inserting the last 10 mm, use a slight twisting and sliding movement to install the annealing tube in the correct position.

The annealing tube must be clean and dry, inside and outside, before it is installed in the sterilization orifice **(4a)**.

The annealing tube is a consumable item. Check it regularly for damage. Replace it if opacity, cracks, chips or similar are observed.  
**Quarz glass annealing tube: Art.-No.: 4.000.400**

Never operate the device without the annealing tube **(4)** in place!

## 6.2 Cleaning the glass front

The glass front can be cleaned with commercial glass cleaners and a soft, lint-free cloth.

Never clean the glass front with abrasive sponges or abrasive cleaners because this may cause damage to the surface.

## 6.3 Cleaning the device

Proceed as in **section 6**. Additional remove power cord from the electrical outlet.

It can be cleaned with a damp cloth and commercially available cleaning agents. Care should be taken when cleaning that no liquid gets inside the unit. It is possible to clean the appliance with a spray-and-wipe disinfectant. The instrument itself must not be sterilized.

Flaming the device is not possible.

Depending on the residual heat, after switching off, the fan **(8)** runs until the excess heat has been dissipated.

## 7. Fan

Smart fan control ensures that even in continuous operation the temperature of the housing does not rise above 65 degrees. As an additional safety mechanism, there is a mechanical thermal circuit breaker.

Do not cover the ventilation grilles **(7)** and cooling fan **(8)** on the unit.

Depending on the residual heat, after switching off, the fan **(8)** runs until the excess heat has been dissipated.

## 8. Error display

This section describes the error display and gives you some advices for solving them.

### 8.1 Over-temperature protection

If the device overheats during the sterilization process, the light will switch off automatically and the stop symbol **(1b)** will flash once. The fan continues to run for some time until the device has cooled down (approx. 10 minutes). After that, switch off the device at the power switch **(13)**. Subsequently the device can be switched on again via the power switch **(13)**.

### 8.2 IR-halogen lamp fault

To determine if the IR-halogen lamp is defective, allow the unit to cool down (approx. 10 minutes) then restart the sterilization process. If the IR-halogen spotlight does not light up within the sterilization time of 10 seconds, it is defective and must be replaced. Switch off the device using the power switch **(13)** on the rear panel.

Before opening the device, remove the power plug **(15)** from the power outlet.

Depending on the residual heat, after switching off, the fan **(8)** runs until the excess heat has been dissipated.

Now the IR-halogen lamp can be changed. The removal and installation instructions are included with the replacement IR-halogen lamp:

**IR-halogen lamp, 230V AC:**

**Art.-No.: 4.001.100**

**IR-halogen lamp, 120V AC:**

**Art.-No.: 4.002.100**

**IR-halogen lamp, 100V AC:**

**Art.-No.: 4.003.100**

### 8.3 Temperature sensor fault

If the temperature sensor has a defect, the sterilization process is stopped and the stop symbol **(1b)** flashes 3 times. Please contact the manufacturer to solve the problem.

## 9. Warranty

All SteriMax basic are covered under our two-year manufacturer warranty against any manufacture defects in material and workmanship (including IR-halogen lamp). Except parts of wear and tear (quartz glass annealing tube). The WLD-TEC warranty guarantees all SteriMax basic under normal usage conditions and does not cover any damages as a direct result of user misuse or / and abuse. The warranty is void upon any unauthorized servicing, disassembly or modifications.

## Troubleshooting guide

EN

- **The Start Symbol (1a) do not light up**

Switch on the device at the power switch **(13)**.

Check if the power plug is correctly wired and make sure that it is connected to a live wall socket (fit to the device).

Check the fuse; if the fuse is faulty, it can be replaced. To do this, pull out the fuse holder **(11)**. Replace the blown fuse with a new one and push the fuse holder back into its socket.

The device is fitted with a (230V: 3,15 AmT; 100V / 120V: 6,30 AmT) protection fuse. The factory provides a spare fuse in the reserve clip in the fuse holder **(11)**.

To replace a fuse, unplug the power cable **(16)** from the power socket **(12)**.

Depending on the residual heat, after switching off, the fan **(8)** runs until the excess heat has been dissipated.

---

- **Inoculating loops are not sterilized**

Annealing tube is dirty. Please clean and dry it, or if necessary replace it **(see section 6.1)**.

---

- **The unit gets hot, the red light turns off during sterilization**

The thermal circuit breaker is switching off the unit **(see section 8.1)**.

Check the functioning of the fan **(see section 7)**.

### Service address:

WLD-TEC GmbH  
Production & Service  
Halle-Kasseler-Str.49  
D-37318 Arenshausen  
Germany

Phone: +49 36081 68940  
Fax: +49 36081 68942  
Email: sales@wld-tec.com  
Internet: www.wld-tec.com

## Technical Data

Touch-Operation:	start of the sterilization (10 sec.) / stop of the sterilization
Foot pedal-Operation*:	start of the sterilization (10 sec.) / stop of the sterilization (*foot pedal is optional)
Status display:	optical display of the sterilization progress
Thermocontrol:	dynamic sterilization timer, reduces the activation of the IR-halogen lamp according to the residual heat
<b>Sterilization at start up:</b>	automatically sterilization after the device is switched on
<b>Annealing tube</b>	
Quartz glass:	750°C - 1000°C (1382°F - 1832°F) sterilization temperature
<b>Safety features</b>	
Residual heat display:	indicates a hot annealing tube
Automatic unit switch off:	after 15 minutes without operation
Temperature monitor:	thermal circuit breaker and fan overrun when switching off the device
<b>Electrical</b>	
<b>Power supply:</b>	<b>Art. No.: 4.101.000:</b> 220 - 240 V 50/60 Hz 2.6 A <b>Art. No.: 4.102.000:</b> 112 - 128 V 50/60 Hz 5.1 A <b>Art. No.: 4.103.000:</b> 95 - 105 V 50/60 Hz 6 A
Power consumption:	approx. 600 W / short time (5 - 10 seconds)
Protection class:	IP 20
Heat source:	IR-halogen lamp
<b>Mechanical</b>	
Quartz glass annealing tube:	outer-Ø 19 mm, inner-Ø 16 mm, length 112 mm
Casing:	heat resistant glass / stainless steel UV and solvent resistant
Measurements (w x h x d):	110 x 170 x 180 mm
Weight:	approx. 1800 g
<b>Information on radiation</b>	
Emission limit:	The emission limit "risk free" is respected. The device is assigned to the risk free group according to EN 62471
<b>Licences</b>	
CE:	EN 61326-1, EN 61010-1, EN 61010-2-010, EN 62471
EU Guidelines:	2014/30/EU, 2014/35/EU, 2011/65/EU, 2006/25/EG
<b>Protected Design:</b>	Ⓢ DE 20 2016 101 067

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**Geschütztes Design:** © DE 20 2016 101 067